



Dongguan Xintai Instrument Co.,Ltd.

PROFESSIONAL INSTRUMENT MANUFACTURER

Body Infrared Thermometer HT-820D

- ◆ Non-contact measuring temperature accurately.
- ◆ It can reduce the risk of virus spread effectively.
- ◆ reading temperature display quickly.response time is only 0.5s.
- ◆ Accuracy is $\pm 0.3^{\circ}\text{C}$ ($\pm 0.6^{\circ}\text{F}$)
- ◆ 34 measuring datas can be saved.



Normal



Slight



Serious

Specification

Storage temperature	-10~60°C(14~140°F)
Usage temperature	10~40°C(50~104°F)
Usage battery	3V
Weight	122g
Size	145 x 80 x 40mm

Technical Parameters

Temperature measuring range	32~43°C(90~109°F)
temperature measuring accuracy	$\pm 0.3^{\circ}\text{C}$ ($\pm 0.6^{\circ}\text{F}$)
Surface Temperature measuring range	0~100°C(32~212°F)
Surface temperature measuring accuracy	$\pm 1^{\circ}\text{C}$ ($\pm 1.8^{\circ}\text{F}$)
Minimum readings	$\pm 0.1^{\circ}\text{C}$ ($\pm 0.1^{\circ}\text{F}$)
Response time	0.5s
Measuring distance	5~8 cm
Auto power off	7s

Different temperature range different backlight, distinguished quickly.



MADE IN CHINA

